

Application Note


Efficient Discovery of B Antigen-Cleaving Enzymes Using Ultrahigh-Throughput Droplet Microfluidics

Charlotte Olagnon, Jacob F. Wardman, and Stephen G. Withers, *et al.* Department of Chemistry, University of British Columbia

Ultrahigh-Throughput Single Emulsion Droplet Screening for the Discovery of New B Antigen Cleaving Enzymes

ACS Catal. 2024, 14, 12884–12894

<https://pubs.acs.org/doi/10.1021/acscatal.4c02165>



Graphical Abstract: Enzyme Screening Strategy

Introduction

ABO blood type mismatch is a critical challenge in blood transfusion. B-type red blood cells (RBCs) carry B antigens, and transfusions from a B-type donor to an incompatible recipient can lead to immune reactions, potentially resulting in death. The goal is to convert B-type RBCs to universal O-type RBCs by removing the terminal sugar (galactose) from the B antigen. Previous enzyme-based methods for this conversion have either been inefficient or required impractical amounts of enzyme. This study presents a **droplet-based microfluidic screening system** to discover new enzymes from metagenomic libraries capable of cleaving the B antigen efficiently.

Methods

A high-throughput single emulsion droplet microfluidic system was developed to encapsulate single cells from metagenomic libraries along with a fluorogenic B antigen substrate within water-in-oil droplets (Fig. 1). Fluorescent droplets, containing cells capable of cleaving the B antigen, are then detected using fluorescence-activated droplet sorting via the On-chip Sort, a microfluidic-based cell sorter. The positive droplets then undergo validation, cloning, and characterization of the discovered enzymes.




Figure 1. Workflow of the screening process: (A) Droplet Generation: Single cells from a metagenomic library are encapsulated in pL-sized droplets along with the B antigen substrate and coupling enzymes. (B) Sorting: Fluorescent droplets are sorted using On-chip Sort (On-chip Biotechnologies) (C) Cloning and Validation: Positive hits are cloned, sequenced, and characterized to confirm their ability to cleave the B antigen.

Results

A fosmid library was created using *E. coli* with metagenomic DNA sourced from the feces of AB blood type donors. This library was encapsulated in droplets and cultured for 48 hours to detect any α -galactosidase activity that cleaves the terminal sugar from the B antigen. From the screening of 1.4×10^5 droplets, 93 (0.07%) positive droplets with a fluorescence exceeding a set threshold (>3-fold background) were collected. After secondary screening of the recovered clones, *PvGH110*, a glycoside hydrolase family 110 α -1,3-galactosidase, was successfully identified. This enzyme was derived from *Phocaeicola vulgatus*, a bacterium in the human gut microbiome.

Subsequent cloning, expression, and purification of the *PvGH110* in *E. coli* allowed for the evaluation of its enzyme activity. When compared to AgaBb and *BfGal110A*, which are known for having the highest specific activities against the B antigen, *PvGH110* demonstrated similar conversion efficiency. Furthermore, in the presence of dextran, which enhances enzyme reactions through macromolecular crowding effects, *PvGH110* was able to completely convert B-type red blood cells into O-type red blood cells (Fig. 2).

The use of droplet technology in this study enhanced the screening efficiency compared to traditional well-plate-based methods and significantly reduced the quantities of enzyme and substrate needed. For instance, only about 30 μ g of substrate was used in this screening, whereas conventional methods would have require approximately 75 mg, marking a 2500-fold reduction. The study presents a robust, efficient, and accessible system for discovering new enzymes capable of converting B-type RBCs to O-type. The identified enzyme, *PvGH110*, is a promising candidate for clinical applications, and the method could be further utilized to discover other antigen-cleaving enzymes, paving the way for future advancements in blood compatibility and transfusion medicine.




Figure 2. Evaluation of B antigen-cleaving enzyme activity: This figure shows the activity of the newly discovered enzyme, *PvGH110*, in converting B⁺ RBCs to O-type RBCs. It compares the conversion efficiency of *PvGH110* with known enzymes (AgaBb, *BfGal110A*), evaluated in the presence and absence of dextran.