

A High-Throughput Droplet-Based Microfluidics Platform for Novel Enzyme Discovery and Engineering

SUMMARY:

In this free webinar, learn about the applications of a droplet-based microfluidics sorter for ultrahigh-throughput enzyme discovery. Attendees will learn the details behind how the microfluidic chip-based cell sorter works. The featured speakers will discuss how the droplet-based sorter can be applied to enzyme discovery and engineering workflows.

BODY:

Droplet-based microfluidics has proven itself to be a powerful tool for the discovery and engineering of enzymes. By miniaturizing assays within droplets to the pL scale, microfluidics enables the screening of large genetic libraries with substantially less substrate use compared to conventional microtitre plates. Moreover, with the appropriate droplet sorter this can be conducted high throughput with libraries of 10⁵ clones being readily screened. In this webinar, the featured speakers will describe how they have applied the On-chip Sort for ultrahigh-throughput enzyme discovery (via functional metagenomic screening) and enzyme engineering by directed evolution.

Both of these techniques allow for the discovery of novel enzymatic activities through exploration of sequence space. In the functional metagenomic screen, they were able to isolate a single active clone from a library of 19,500 clones in a single sort. Such low hit rates are typical in functional metagenomics and highlight the need for high-throughput assays in order to identify active clones. Meanwhile, a single round of the directed evolution campaign in droplets (with a library of 100,000 clones) revealed a panel of clones with improvements from five to 35-fold enhancements in activity compared to the parent.

Both of these examples are enabled by the miniaturization offered by microfluidics as the screens required the use of oligosaccharide substrates. At these scales the cost of substrate would be prohibitive if done in microtitre plates. Such results highlight the utility of droplet-based microfluidics for the rapid exploration of sequence space where other assays may be otherwise unfeasible.

Join this webinar to learn about the applications of a droplet-based microfluidics sorter for ultrahigh-throughput enzyme discovery.

Join experts from the University of British Columbia, Jacob Wardman, PhD Candidate; and Charlotte Olagnon, Postdoctoral Research Fellow, for the live webinar on Wednesday, March 1, 2023, at 1pm EST (10am PST).

For more information, or to register for this event, visit [A High-Throughput Droplet-Based Microfluidics Platform for Novel Enzyme Discovery and Engineering](#).